

FIG. - 1




FIG. - 2





FIG.-4A




FIG.-4B




FIG. - 7





FIG. - 5B




FIG. - 5C

FIG. - 6





FIG. - 8A

FIG. - 8B

FIG. - 8C

208		200		202		204	
BC	Selected bc (6.9-10.4/0.1) (7.70-9.1/.05)	8.40	Suggested Base Curve is 8.4				
J1	Radial distance (OZ2) from the lens center to 1st junction mm (1.0-5.9/0.1)	210	3.00	1	corneal apical radius (mm)	7.58	
SW	Width of the S curve mm (.75,1)	1.00	EYE				
MAT	Lens material (FP30, FP60, FP92, FPI51, HDS, Other)	212		Ref. Index of material used = 1.449 If 'other' was selected input RI in Cell H4			
P	lens power desired (-1.00, -0.50, 0.00, 0.50, 0.75, 1.0, 1.25, 1.5, 2.0, 2.5)	214		Front Surface central radius = 8.37			
92	Delta R (mm) translation of 1st junction radially from BC origin (0.08-0.2/0.02)	222	0.14	True center thickness (mm) = 0.152			
92	Delta R (mm) translation of 2nd junction radially from BC origin (0.1-0.22/0.02)	0.18		true offset between landing zones at J2 = 0.179			
A	Angle of the landing zone (-25.5 to -50.0/.5)	216	-35.00	Present lens height (mm) above cornea at diameter of tangential touch = 0.040			
D	selected lens diameter mm (8.0-12.9/0.1)	209	10.50	Diameter recommended from HVID = 10.6			
SD	Selected depth of the S curve mm (.15-1.0/.05) (0.3-0.65/.025) use next smaller than est.	0.500		Recommended depth (mm) S curve for desired correction @6u/D = 0.510 mm			
		220	218				
				lens / cornea power (D) difference wanted			
				Actual power (D) difference between bc and apical cornea = -4.35			
				Recommended diameter for lentic = 8.024			
				recommended radius of curve for lentic = 8.106			
				Origin for lentic curve is on y axis displaced from apex of front curve = 8.068			
				Estimated elevation at J2 = 0.070			
				Diameter where LZ would make tangential touch = 9.08			
				Dia giving desired LZ lift = 10.42			
				Edge lift at selected diameter = 0.094			
				ellipticity of the cornea			
				Desired edge lift (mm) when landed at full Diameter = 0.083			
				Ab, the long axis of the ellipse creating the base curve edge (below)			
				FOR SPHERICAL FRONTS target edge thickness below			
				0.40			
				Af, the long axis of the ellipse creating the front curve edge (below)			
				SPHERICAL FRONTS- max thickness peripheral to J1 before lentic (in mm > Delta 2) see below			
				0.40			
				base to front at which the transition from base ellipse to front ellipse is found (below)			
				Minimum thickness peripheral to J1 before lentic (in mm > Delta 1) see below			
				0.25			
				0.006			
				0.01			

FIG.-10




FIG. -11




FIG.-12





FIG.-13

FIG-14

FIG-14

FIG.-15




FIG.-16




FIG.-17





FIG.-18



Selected bc		8.35		Suggested Base Curve is 8.3		200		202		204	
BC	(6.9-10.4/0.1) (7.70-9.1/0.05)	2.50		3B		210		202		204	
J1	Radial distance (OZ2) from the lens center to 1st junction mm (1.0-5.9/0.1)	2.50		3B		210		202		204	
SW	Width of the S curve mm (.75,1)	2.00		EYE		212		202		204	
MAT	Lens material (FP30, FP60, FP92, FP151, HDS, Other)	HDS		Ref. Index of material used = 1.449 If 'other' was selected input RI in Cell H4		7.25		lens / cornea power (D) difference wanted		ellipticity of the cornea	
P	lens power desired (-1.00, -0.50, 0.00, 0.50, 0.75, 1.0, 1.25, 1.5, 2.0, 2.5)	214		Front Surface central radius = 8.32		Volume between BC and cornea (uL) = 0.699		Actual power (D) difference between bc and apical cornea = -6.13		Desired edge lift (mm) when landed at full Diameter = 0.09272	
Q2	Delta R (mm) translation of 1st junction radially from BC origin (0.08-0.2/0.02)	0.14		True center thickness (mm) = 0.148		Volume between S curve and cornea (uL) = 2.812		Recommended diameter for lentic = 6.784		Ab, the long axis of the ellipse creating the base curve edge (below)	
A	Angle of the landing zone (-25.5 to -50.0/5)	-38.00		Present lens height (mm) above cornea at diameter of tangential touch = 0.038		Volume between pretouch Landing Zone and cornea (uL) = 0.122		recommended radius of curve for lentic = 7.615		FOR SPHERICAL FRONTS target edge thickness below	
D	selected lens diameter mm (8.0-12.9/0.1)	10.20		Diameter recommended from HVID = 10.2		TOTAL VOLUME = 3.633(uL)		Origin for lentic curve is on y axis displaced from apex of front curve = 7.541		SPHERICAL FRONTS- max thickness peripheral to J1 before lentic (in mm > Delta 2) see below	
SD	Selected depth of the S curve mm (.15-1.0/0.05) (0.3-0.65/ .025) use next smaller than est.	1.116		Recommended depth (mm) S curve for desired correction @6u/D = 1.116 mm		Diameter where LZ would make tangential touch = 9.21		Estimated elevation at J2 = 0.040		Minimum thickness peripheral to J1 before lentic (in mm > Delta 1) see below	
						Dia giving desired LZ lift = 10.53		fixed (rear thickness)			
						Edge lift at selected diameter = 0.071		0.006		0.25	
										0.01	

FIG.-20




FIG.-21




FIG.-22





FIG.-23



[illegible]

FIG-24

FIG.-25




FIG.-26




FIG.-27




FIG.-28





FIG.-29



Selected bc		Suggested Base Curve is 8.4					
BC	(6.9-10.4/0.1) (7.70-9.1/0.05)	8.40					
J1	Radial distance (OZ2) from the lens center to 1st junction mm (1.0-5.9/0.1)	3.00	5B		corneal apical radius (mm)	lens / cornea power (D) difference wanted	ellipticity of the cornea
SW	Width of the S curve mm (.75,1)	1.00	EYE		7.75	-3.50	0.7
MAT	Lens material (FP30, FP60, FP92, FP151, HDS, Other)	HDS		Ref. Index of material used = 1.449 If 'other' was selected input RI in Cell H4		Actual power (D) difference between bc and apical cornea = -3.37	Desired edge lift (mm) when landed at full Diameter = 0.077
P	Lens power desired (-1.00, -0.50, 0.00, 0.50, 0.75, 1.0, 1.25, 1.5, 2.0, 2.5)	0.50		Front Surface central radius = 8.36	Volume between S curve and cornea (uL) = 1.195	Recommended diameter for lentic = 7.735	Ab, the long axis of the ellipse creating the base curve edge (below)
Q1	Delta R (mm) translation of 1st junction radially from BC origin (0.08-0.2/0.02)	222 0.10		True center thickness (mm) = 0.110	Volume between pretouch Landing Zone and cornea (uL) = 0.439	recommended radius of curve for lentic = 9.295	
Q2	Delta R (mm) translation of 2nd junction radially from BC origin (0.1-0.22/0.02)	242 0.10		true offset between landing zones at J2 = 0.100	TOTAL VOLUME = 2.382(uL)	Origin for lentic curve is on y axis displaced from apex of front curve = 9.400	SPHERICAL FRONTS- max thickness peripheral to J1 before lentic (in mm > Delta 2) see below
A	Angle of the landing zone (-25.5 to -50.0/.5)	-32.50		Present lens height (mm) above cornea at diameter of tangential touch = 0.027	Diameter where LZ would make tangential touch = 8.99	Estimated elevation at J2 = 0.047	
D	selected lens diameter mm (8.0-12.9/0.1)	10.00		Diameter recommended from HVID = 10	Dia giving desired LZ lift = 10.59	fixed (tear thickness)	base to front at which the transition from base ellipse to front ellipse is found (below)
SD	Selected depth of the S curve mm (.15-1.0/.05) (0.3-0.65/.025) use next smaller than est.	0.475		Recommended depth (mm) S curve for desired correction @6w/D = 0.478 mm	Edge lift at selected diameter = 0.048	0.006	Minimum thickness peripheral to J1 before lentic (in mm > Delta 1) see below
						0.25	0.01

FIG-30




FIG.-32




FIG.-33





FIG.-34



Suggested Base Curve is 7.5													
BC	Selected bc (6.9-10.4/0.1) (7.70-9.1/0.05)	7.50											
J1	Radial distance (OZ/2) from the lens center to 1st junction mm (1.0-5.9/0.1)	2.50	5B			corneal apical radius (mm)	7.8	lens / cornea power (D) difference wanted	2.00	ellipticity of the cornea	0.3	HVID (mm)	11.9
SW	Width of the S curve mm (.75,1)	1.50	EYE										
MAT	Lens material (FP30, FP60, FP92, FP151, HDS, Other)	2.10	HDS	Ref. Index of material used = 1.449 If "other" was selected input RI in Cell H4		Volume between BC and cornea (uL) = 0.298		Actual power (D) difference between bc and apical cornea = 1.73		Desired edge lift (mm) when landed at full Diameter = 0.062			1.45
P	lens power desired (-1.00, -0.50, 0.00, 0.50, 0.75, 1.0, 1.25, 1.5, 2.0, 2.5)	0.50		Front Surface central radius = 7.49		Volume between S curve and cornea (uL) = 1.383		Recommended diameter for lentic = 5.737		Ab, the long axis of the ellipse creating the base curve edge (below)		FOR SPHERICAL FRONTS target edge thickness below	
Q1	Delta R (mm) translation of 1st junction radially from BC origin (0.08-0.2/0.02)	0.14		True center thickness (mm) = 0.149		Volume between pretouch Landing Zone and cornea (uL) = 0.491	246	recommended radius of curve for lentic = 8.482			0.40		0.18
Q2	Delta R (mm) translation of 2nd junction radially from BC origin (0.1-0.22/0.02)	0.18		true offset between landing zones at J2 = 0.180		TOTAL VOLUME = 2.171(uL)		Origin for lentic curve is on y axis displaced from apex of front curve = 8.553		Af, the long axis of the ellipse creating the front curve edge (below)		SPHERICAL FRONTS- max thickness peripheral to J1 before lentic (in mm > Delta 2) see below	
A	Angle of the landing zone (-25.5 to -50.0/5)	-35.00		Present lens height (mm) above cornea at diameter of tangential touch = 0.024		Diameter where LZ would make tangential touch = 9.08		Estimated elevation at J2 = 0.056			0.40		0.01
D	selected lens diameter mm (8.0-12.9/0.1)	10.00		Diameter recommended from HVID = 10.9		Dia giving desired LZ lift = 10.52		fixed (tear thickness)		base to front at which the transition from base ellipse to front ellipse is found (below)		Minimum thickness peripheral to J1 before lentic (in mm > Delta 1) see below	
SD	Selected depth of the S curve mm (.15-1.0/0.05) (0.3-0.65/.025) use next smaller than est.	0.636		Recommended depth (mm) S curve for desired correction @6w/D = 0.646 mm		Edge lift at selected diameter = 0.047		0.024			0.25		0.01

FIG-35

FIG.-36




FIG. 37




FIG.-38




FIG.-39




FIG. - 40

Determine necessary spectacle correction, add up to 0.5 diopter overcorrection to extend duration of treatment

350

Determined corneal curvature from keratometry or topography and add 0.2 mm per Diopter of needed correction to select lens base curve.

352

Measure horizontal visible iris diameter and select a lens diameter approximately 1 mm smaller than HVID

354

Use topography information or trial fitting with various landing zones, to find the angle whose point of tangency is half way between the selected diameter and 8mm (sum of standard central zone width (6mm) and connector zone width 1mm (2mm considering both sides).

356

Using fluoresce drops with trial lenses (of the selected diameter and base curve and angle and having various connector zone depths) estimate the rectangle depth that leaves the tangent point elevated above the cornea 6 microns per dioter of needed correction. Similar estimates can be made using topographical information.

358

Compare to model eye using software and confirm lens fit on patient. Review LZ elevation, edge lift, apical contact and centration.

360

458

FIG.-4

FIG-42

